

Combined Pharmacy Refrigerator and Freezer



Suitable for hospital, blood-bank, center for Disease Prevention, clinic and scientific research institution, etc.

Used for refrigeration preserve vaccines, reagents, etc.

Characteristics of the product function:

■Features Refrigeration System

Optimized refrigeration technology, more powerful refrigerating capacity;
CFC-Free foaming, CFC-Free refrigerant, more environmentally friendly;

■Temperature Control

Equipped with micro-computer controller, Refrigerator and freezer temperature can be display separately at the same time;

The temperature in the refrigerating chamber can be adjusted from:2C to 8C,
freezing chamber temperature: -15--26C;

It adopts forced ventilation refrigeration, which can ensure the temperature in the refrigerating chamber is steady.

■ Safety Control

Multiple alarming methods high temperature alarm, low temperature alarm, sensor malfunction alarm, open the door alarm;

Two alarming methods: buzzer alarming, light flashing alarming; Protection function: power-on delay, time between stops;

■ Ergonomic Design

LCD liquid crystal display, which is convenient to observe;

Refrigerator chamber equipped with plastic impregnated steel wire shelves and freezer chamber with transparent drawers.

Safety door lock preventing unauthorized access;

Wide voltage design from 192V to 242V;

Technical Parameters

Model		YJ-PCF265	YJ-PCF288
Technical Data	Cabinet Type	Vertical	Vertical
	Climate Class	ST	ST
	Cooling Type	Forced air cooling	Forced air cooling
	Defrost Mode	Auto	Auto
	Refrigerant	CFC-Free	CFC-Free
Performance	cooling performance(℃)	R: 4/F:-25	R: 4/F:-40
	Temperature Range(℃)	R: 2~8/F:-10~-25	R: 2~8/F:-10~-40
Control	Controller	Microprocessor	Microprocessor
	Display	LED	LED
Material	Interior	ABS	ABS
	Exterior	Galvanized steel powder coating	Galvanized steel powder coating
	Door	Foaming door	Foaming door

Electrical Data	Power Supply (V/Hz)	220/50	220/50
	Power (W)	90	290
	Electrical Current (A)		
Dimensions	Capacity (L)	265	265
	Net/Gross Weight (approx)	60/65	60/65
	Interior Dimensions (W*D*H)	R: 460x480x720mm; F: 405x460x692mm	R: 588x552x630mm; F: 511x491x415mm
	Exterior Dimensions (W*D*H)	575x641x1687mm	737x780x1860mm
	Packing Dimensions (W*D*H)	660x700x1723mm	850x800x1950mm
	Container load (20/40' /40' H)	27/54/54	14/30/30
Functions	High/Low Temperature	Y	Y
	Sensor Failure	Y	Y
	Door Ajar	Y	Y
	Lockage	Y	Y
	Inner LED light	Y	Y
Accessories	Caster	Y	Y
	Shelves/Drawer	3/4	3/3